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## General Information

These instructions for use must be read carefully every time before the products are used. They must be provided to anyone entrusted with the care of patients for informational purposes.

### 1. Intended use

Thoracic and trocar catheters are employed in cases where air and fluid need to be drained out of the pleural cavity and/or mediastinum under medical circumstances.

### 2. Product description

Thoracic catheters and trocar catheters are used in conjunction with a thoracic drainage system to collect and remove fluids and/or air from the mediastinal or pleural cavity following injury- or surgery-related opening of the thorax. Every catheter is equipped with a radiopaque stripe to aid radiological visual control. In order to facilitate positioning and removal of the catheter, depth markings are included at regular instances with the scale starting from the last proximal catheter eye.

#### Thoracic catheter (straight/curved)

The curved version of the thoracic catheter is designed in such a way as to be positioned between the visceral pleura and parietal pleura, with the 90° angle allowing proximal or distal positioning. A special, conical connector has been designed to allow the thoracic catheter to be pulled through the thoracic wall atraumatically. The connector is cut following positioning of the catheter so as to ensure secure attachment to the drainage tube. The extra-long thoracic catheter is often used in combination with standard length thoracic catheters in order to ensure locally differentiated (or multilocular) drainage.

#### Trocar catheter

The conical puncture tip has been designed in such a way as to allow the catheter to be positioned through the thoracic wall as atraumatically as possible. The distal end of the catheter tube is equipped with a connection adapter to which suitable drainage systems can be connected. There is a metal trocar in the catheter tube extending into the puncture tip for improved guiding. The distal end of the trocar features a knob-shaped handle.

### 3. Period of use

Replace the thoracic or trocar catheter after 7 days at the latest or sooner if it becomes damaged or blocked.

### 4. Indication

The catheter is intended for short-term use in cases where the medical circumstances require the drainage of air and/or liquid from the pleural cavity and/or mediastinum.

Thoracic and trocar catheters are used

- in patients in whom symptomatic conditions have resulted in air entering the pleural cavity
- in patients with bleeding, fluid accumulation and purulent effusions within the thorax, and
- in patients following surgical interventions in the thoracic cavity or following heart surgery.

## 5. Contraindications

There are no known contraindications.

## 6. Side effects

There are no known side effects.

## 7. Warnings /Safety precautions

- Prior to positioning of the thoracic/trocar catheter, the introduction site is located on the thoracic wall and a catheter of the corresponding size selected.
- The attending doctor is responsible for selecting the size of the thoracic/trocar catheter to be used based on what is expected to be drained.
- Generally, catheters for draining air are positioned anteriorly and high whilst catheters for draining fluids are positioned posteriorly and low.
- Ensure that the thoracic/trocar catheter is always patent and not kinked.
- The ingress of air must be avoided when removing the catheter. An air-tight dressing should thus be applied immediately (purse-string suture if applicable).
- The thoracic/trocar catheter can be withdrawn for correction but not subsequently advanced so as to avoid infection. A new thoracic/trocar catheter should not be inserted through the same access for hygienic reasons.

## 8. Use

### 8.1. Positioning

#### Thoracic catheter

1. The sterile thoracic catheter is removed from its protective packaging maintaining aseptic conditions.
2. A stab incision is made at the exit site on the skin.
3. The catheter can then be positioned in the pleural cavity through the incision from the outside in or from the inside out. The doctor makes the decision depending on the indications in each case.
4. Following correct positioning, the blue cap (not present in silicone versions) on the end of the catheter is cut (the piece cut off must be disposed of) and a prepared and suitable drainage system connected.
5. Securing of the connection between the connector and drainage tube and the catheter itself with medical tape is recommended.
6. The catheter should be fixed to the skin with a suture.

#### Trocar catheter

1. The sterile trocar catheter is removed from its protective packaging maintaining aseptic conditions.
2. Once the incision has been made, the pleural cavity is punctured with the trocar catheter and the catheter advanced to the desired position.
3. Once correctly positioned, the trocar is removed from the catheter tube.
4. The catheter is then connected to a prepared and suitable drainage system.
5. Securing of the connection between the connector and drainage tube and the catheter itself with medical tape is recommended.
6. The catheter should be fixed to the skin with a suture.

### 8.2. Removal

1. Any fixation points (purse-string suture) and medical tape should be removed and then

the catheter closed off with a clamp close to the patient to prevent air from entering the thoracic cavity.

2. Remove the drainage system.
3. Pull out the catheter. When doing so, take care to avoid ingress of air.
4. An air-tight dressing should be applied immediately.

## **9. Cleaning/ Disinfection**

No cleaning or disinfection.

## **10. Storage/ Warehousing**

The catheters should always be stored in a dry environment protected from dust, direct exposure to sunlight and/or heat.

## **11. Disposal**

Disposal must be performed in accordance with the applicable national regulations for waste materials.

## **12. Reporting**

All serious incidents and events occurring in connection with the product must be reported to the manufacturer and the responsible body.

## **13. Guarantee**

We guarantee that the delivered object is free from defects in accordance with the scope of our general terms and conditions (GTCs). You can find a copy of our GTCs among our business documentation and on our website at [www.primed-halberstadt.de](http://www.primed-halberstadt.de), etc.

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